F5

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EN 1. Warning
HEALTH RISK
Can cause irritation during contact with skin and eyes. Prolonged exposure to the product can cause allergic reactions.

USE PRECAUTIONS
Do not use if the product or the packaging is damaged, or if the product has changed in appearance. Do not use after expiration date. Comply with the proportions of 40,0 gr of F5 per litre, otherwise there is a risk of reducing the performance of the preservation medium. Do not store close to food, beverages or animal feed. Use gloves and a protective mask when handling the product.

STORAGE
The medium can be transported and preserved at room temperature not exceeding +25°C. Do not freeze, either in powder or in solution form.

2. Use protocol
PREPARATION
To prepare the final solution it is necessary to dilute 40 gr of F5 for every liter of apyrogenic sterile water or type II.

It is important to measure the volume of sterile pyrogen-free water accurately to avoid any variation in osmotic pressure.

USE
- Prepare the extender at least half an hour prior to use for equilibration.

- Verify the temperature of the semen and that of the preservation medium. Do not mix the preservation medium and the semen if the difference in temperature between the two is more than 2° C.

- Use reconstituted medium within 48 hours maximum.
- Any unused portion of the reconstituted medium after 6 hours should be stored in a refrigerator between +2° and +8°C in a close container.

- The extender cannot be frozen, in either powder or solution form.
DILUTION RECOMMENDATION
We recommend making a 1:1 pre-dilution (interval between collection - pre-dilution less than 15 minutes) with extender kept between +32°C and +34°C.

- Dilute to the final proportions to obtain 2.5 and 3 billion living spermatozooids per dose of 100 ml depending on the original quality of the ejaculate and the storage conditions of the semen.

- In the absence of a count or photometric measurement of the spermatozoid concentration, the minimum dilution to be respected is 1 volume of semen for 4 volumes of extender with a maximum of 12 doses of 100 ml.

- Store the doses between +15°C and +17°C in the dark.
- Use the doses within 5 days of dilution.

3. Evaluation of samples of diluted semen
Preparation of samples
Place 10 ml of diluted semen taken from a dose preserved at 17°C in a culture tube with a leak-tight seal.

Seal the tube.
Mix gently the tube by turning it over 3 times in succession.
Incubate for 10 minutes at 37°C in a water bath.

EVALUATION OF SAMPLES
We recommend using pre-calibrated microscope slides (Leja type). Otherwise, use new clean microscope slides and cover slides.

Warm the slides at 37°C using a slide warmer.
Open the tube previously warmed up to 37°C, and place the semen in a pre-calibrated slide chamber.

Make an immediate assessment of the quality of the semen.
4. Standard of production and quality control
Physical-chemical characteristics of the reconstituted medium (with water)

- Osmolarity: 290 +/- 20 mOsm/kg H₂O
- pH: 7.1 +/- 0.2

PRODUCTION STANDARD
The extender is manufactured in a MEDI NOVA laboratory, under hygrometry and temperature controlled atmosphere. MEDI NOVA is certified ISO 9001 2008. The extender is tested by the quality department of MEDI NOVA for compliance with established specifications.

CERTIFICATE
A batch certificate is available upon request.

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